


## Disposable Grasping Forceps



Disposable Grasping Forceps  
【Rat Tooth Type】



Disposable Grasping Forceps  
【Rat Tooth with Alligator Type】



Disposable Grasping Forceps  
【Pelican Type】



Disposable Grasping Forceps  
【Alligator Type】



Disposable Grasping Forceps  
【3-Prong Type】



Disposable Grasping Forceps  
【5-Prong Type】



Disposable Grasping Forceps  
【2-Prong Type】

### Disposable Grasping Forceps 【Rat Tooth Type】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810GI	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1812GI	Ø1.8x1200mm	≥Ø2.0mm	Bronchoscopy
AF-D1816GI	Ø1.8x1600mm	≥Ø2.0mm	Gastroscopy
AF-D2416GI	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2418GI	Ø2.4x1800mm	≥Ø2.8mm	Gastroscopy
AF-D2421GI	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GI	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy

Rat tooth for flat and hard objects ( e.g.coins ).

Opening width:  
5mm(Ø1.8) 6mm(Ø2.4)



### Disposable Grasping Forceps 【Rat Tooth with Alligator Type】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810GJ	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1812GJ	Ø1.8x1200mm	≥Ø2.0mm	Bronchoscopy
AF-D1816GJ	Ø1.8x1600mm	≥Ø2.0mm	Gastroscopy
AF-D2416GJ	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2418GJ	Ø2.4x1800mm	≥Ø2.8mm	Gastroscopy
AF-D2421GJ	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GJ	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy

Sharp tooth, long jaws for difficult retrieval of flat and hard objects.

Opening width:  
5mm(Ø1.8) 6mm / 8.5mm(Ø2.4)




### Disposable Grasping Forceps 【Pelican Type】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810GK	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1812GK	Ø1.8x1200mm	≥Ø2.0mm	Bronchoscopy
AF-D1816GK	Ø1.8x1600mm	≥Ø2.0mm	Gastroscopy
AF-D2416GK	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2418GK	Ø2.4x1800mm	≥Ø2.8mm	Gastroscopy
AF-D2421GK	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GK	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy

For big objects. Opening width: 15mm.

Opening width:  
9.5mm(Ø1.8) 15mm(Ø2.4)



### Disposable Grasping Forceps 【Alligator Type】


Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810GH	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1812GH	Ø1.8x1200mm	≥Ø2.0mm	Bronchoscopy
AF-D1816GH	Ø1.8x1600mm	≥Ø2.0mm	Gastroscopy
AF-D2416GH	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2418GH	Ø2.4x1800mm	≥Ø2.8mm	Gastroscopy
AF-D2421GH	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GH	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy

For sharp-edged objects.


Opening width:  
7mm(Ø1.8) 7.5mm(Ø2.4)



Disposable Grasping Forceps [3-Prong Type]			
Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1807GM103	Ø1.8x700mm	≥Ø2.0mm	Choledochoscopy
AF-D1810GM103	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1815GM103	Ø1.8x1500mm	≥Ø2.0mm	Gastroscopy
AF-D2416GM103	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2417GM103	Ø2.4x1700mm	≥Ø2.8mm	Gastroscopy
AF-D2421GM103	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GM103	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy




Disposable Grasping Forceps [3-Prong Hook Type]			
Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1807GM103-H	Ø1.8x700mm	≥Ø2.0mm	Choledochoscopy
AF-D1810GM103-H	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1815GM103-H	Ø1.8x1500mm	≥Ø2.0mm	Gastroscopy
AF-D2416GM103-H	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2417GM103-H	Ø2.4x1700mm	≥Ø2.8mm	Gastroscopy
AF-D2421GM103-H	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GM103-H	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy




Disposable Grasping Forceps [5-Prong Type]			
Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1807GM105	Ø1.8x700mm	≥Ø2.0mm	Choledochoscopy
AF-D1810GM105	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1815GM105	Ø1.8x1500mm	≥Ø2.0mm	Gastroscopy
AF-D2416GM105	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2417GM105	Ø2.4x1700mm	≥Ø2.8mm	Gastroscopy
AF-D2421GM105	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GM105	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy



Disposable Grasping Forceps [5-Prong Hook Type]			
Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1807GM105-H	Ø1.8x700mm	≥Ø2.0mm	Choledochoscopy
AF-D1810GM105-H	Ø1.8x1000mm	≥Ø2.0mm	Bronchoscopy
AF-D1815GM105-H	Ø1.8x1500mm	≥Ø2.0mm	Gastroscopy
AF-D2416GM105-H	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2417GM105-H	Ø2.4x1700mm	≥Ø2.8mm	Gastroscopy
AF-D2421GM105-H	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GM105-H	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy




Disposable Grasping Forceps [2-Prong Type]			
Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D2416GM102	Ø2.4x1600mm	≥Ø2.8mm	Gastroscopy
AF-D2417GM102	Ø2.4x1700mm	≥Ø2.8mm	Gastroscopy
AF-D2421GM102	Ø2.4x2100mm	≥Ø2.8mm	Colonoscopy
AF-D2423GM102	Ø2.4x2300mm	≥Ø2.8mm	Colonoscopy




## Disposable Retrieval Baskets



Disposable Grasping Forceps  
【Basket without Tip Type,  
Nitinol】



Disposable Grasping Forceps  
【Fold Angular Basket Type,  
Stainless Steel】



Disposable Grasping Forceps  
【Spirality Basket Type, Nitinol】