


Disposable Biopsy Forceps



Disposable Biopsy Forceps
【Oval Cups】



Disposable Biopsy Forceps
【Oval Cups With Needle】



Disposable Biopsy Forceps
【Oval Cups Coated】



Disposable Biopsy Forceps
【Oval Cups With Needle Coated】



Disposable Biopsy Forceps
【Alligator Cups】



Disposable Biopsy Forceps
【Alligator Cups Coated】



Disposable Biopsy Forceps
【Alligator Cups With Needle】



Disposable Biopsy Forceps
【Alligator Cups With Needle Coated】

Disposable Biopsy Forceps 【Oval Cups】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BT	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BT	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BT	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BT	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BT	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BT	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BT	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy
AF-D3023BT	Ø3.0x2300mm	≥ Ø3.2mm	Colonoscopy

Oval and fenestrated cups for deeper biopsy.

Opening width:
5mm(Ø1.8) 6.5mm(Ø2.4) 8.5mm(Ø3.0)



Disposable Biopsy Forceps 【Oval Cups Coated】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BTC	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BTC	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BTC	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BTC	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BTC	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BTC	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BTC	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy
AF-D3023BTC	Ø3.0x2300mm	≥ Ø3.2mm	Colonoscopy

Oval and fenestrated cups for deeper biopsy.

Opening width:
5mm(Ø1.8) 6.5mm(Ø2.4) 8.5mm(Ø3.0)



Disposable Biopsy Forceps 【Oval Cups With Needle】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BW	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BW	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BW	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BW	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BW	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BW	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BW	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy

Oval, fenestrated cups with spike for accurate targeting. No sliding for gathering bigger tissue samples.

Opening width:
5mm(Ø1.8) 6.5mm(Ø2.4)



Disposable Biopsy Forceps 【Oval Cups With Needle Coated】

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BWC	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BWC	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BWC	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BWC	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BWC	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BWC	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BWC	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy

Oval, fenestrated cups with spike for accurate targeting. No sliding for gathering bigger tissue samples.

Opening width:
5mm(Ø1.8) 6.5mm(Ø2.4)



Disposable Biopsy Forceps [Alligator Cups]

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BU	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BU	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BU	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BU	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BU	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BU	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BU	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy
AF-D2823BU	Ø2.8x2300mm	≥ Ø3.0mm	Colonoscopy

Serrated, alligator cups for deep biopsy of hard and difficult-to-access tissue.



Disposable Biopsy Forceps [Alligator Cups Coated]

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BUC	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BUC	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BUC	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BUC	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BUC	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BUC	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BUC	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy
AF-D2823BUC	Ø2.8x2300mm	≥ Ø3.0mm	Colonoscopy

Serrated, alligator cups for deep biopsy of hard and difficult-to-access tissue.



Disposable Biopsy Forceps [Alligator Cups With Needle]

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BZ	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BZ	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BZ	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BZ	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BZ	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BZ	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BZ	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy

Serrated, alligator cups for deep biopsy of hard and difficult-to-access tissue. No slippage with spike for accurate targeting.



Disposable Biopsy Forceps [Alligator Cups With Needle Coated]

Item No.	Tube Diameter & Working Length	Working Channel Diameter	Use
AF-D1810BZC	Ø1.8x1000mm	≥ Ø2.0mm	Bronchoscopy
AF-D1812BZC	Ø1.8x1200mm	≥ Ø2.0mm	Bronchoscopy
AF-D1816BZC	Ø1.8x1600mm	≥ Ø2.0mm	Gastroscopy
AF-D2416BZC	Ø2.4x1600mm	≥ Ø2.8mm	Gastroscopy
AF-D2418BZC	Ø2.4x1800mm	≥ Ø2.8mm	Gastroscopy
AF-D2421BZC	Ø2.4x2100mm	≥ Ø2.8mm	Colonoscopy
AF-D2423BZC	Ø2.4x2300mm	≥ Ø2.8mm	Colonoscopy

Serrated, alligator cups for deep biopsy of hard and difficult-to-access tissue. No slippage with spike for accurate targeting.



Grasping Forceps & Retrieval Baskets



DISPOSABLE